

# HydROC Fiber and Hydrolysis Unit


One instrument for both Fiber and Hydrolysis  
Easy to use Fiber bag system  
All samples separated



*specifications subject to change without notice*



*The sample is weighed into an OPSIS LiquidLINE Fiberbag.*



*Sealing of bag with Thermosealer ensures that sample remains inside the bag during the whole process*



*Boiling according to Weende Fiber method.*



*Efficient water rinsing, using the HydROC drain valves.*



*Removal of Filterbag for ashing and weighing.*

OPSIS AB, Box 244  
SE-244 02 Furulund Sweden  
Telephone +46 46 72 25 00  
Telefax +46 46 72 25 01  
E-mail info@opsis.se  
www.liquidline.se



LP1023  
2021 02

**OPSIS**<sup>®</sup>  
LIQUIDLINE<sup>®</sup>

OPSIS AB • +46 46 72 25 00 • info@opsis.se • www.liquidline.se